

Tripod

Documentation

Attaching the light bar to the Tripod

The only accessory you will need for your lightbar is the “Mounting Plate”

First, locate the mounting plate, it will be included in your Simple EMDR Set either in your lightbar carry case OR mounted on the top of the tripod.

The mounting plate allows you to quickly and easily attach or detach the lightbar from the tripod using the quick release lever on the top of the trip



Quick release lever being used to remove mounting plate from tripod.



Lightbar ready to attach mounting plate.



Showing the mounting plate after screwing it to the bottom of the lightbar.



Finally, the lightbar assembly put in place using the quick release lever on the top of the tripod.